Keno™lux S 100

Version: B Date of Update: 02/03/2017



Sanitary cleaner

Product description

Keno™lux S 100 is a concentrated acid cleaning product for daily or periodic cleaning of sanitary spaces.

Properties

Keno™lux S 100 removes and prevents lime scale and fat residues.

Keno™lux S 100 can be used for the removal of fingerprints, soap, residues,...

Keno™lux S 100 is based on lactic acid and removes lime scale and urine deposits.

Keno™lux S 100 has a pleasant perfume and gives a shiny result.

Direction of use

Dilute 0,1-0,2% Keno™lux S 100 in a bucket of water. Rinse after application with clean water.

Spray:

Use a dosage of 1-2%.

Specifications

Physical state Viscous liquid

Colour Pink

рΗ ≈ 2.5 (100%)

Freezing point -3 °C 100 °C **Boiling point** Relative density 1.04 Solubility in water 100 %

Opslag en verwijdering

Keep only in the original container in a cool, well ventilated place. Keep container closed when not in use. Dispose in a safe manner in accordance with local/national regulations.

Safety instructions for the user

Avoid contact between alkaline and acid products.

Always test the product first on a small surface.

Consult your Cid Lines sales' person for a suitable cleaning- and disinfecting program that is adjusted according to the contamination and the specific conditions of the company.

Waterpoortstraat 2, 8900 leper (Belgium) • T +32 (0)57 21 78 77 • F +32 (0)57 21 78 79 • info@cidlines.com • www.cidlines.com

Before use, please check the MSDS.